



TECH BULLETIN

039
03/85

Published regarding engineering changes and improvements

SUBJECT: Flat Valves for Model 390

Beginning with serial number 2850101, the flat valve assembly became standard in the Model 390 pumps.

The Flat Valve parts are not interchangeable with the Quiet Valve parts and must be replaced as a set.

| | OPTIONAL | STANDARD |
|---------------------|---------------|--------------|
| | [Quiet Valve] | [Flat Valve] |
| Valve..... | 43721 | 104302 |
| Seat..... | 43779 | 103956 |
| Retainer | 43780 | 103957 |
| Spring | 43751 | 43751 |
| O-Ring, D.V.S. | 43793 | 43793 |
| Valve Kit..... | 30859 | 34010 |

CAT PUMPS
Technical Services Department

03/96